Retainer

INSTRUCTIONS FOR USE

CONTENTS

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English Version

1. Model and Specification

This manual applies to Retainer. This product is personalized based on dental impressions (physical or digital) of the teeth. The size and shape will vary as a result.

2. Main Structure and Composition of the Product

This product consists of a retainer and an attachment template (optional). The retainer is made of thermoplastic polyurethane (TPU) material, and the attachment template of thermoplastic TPU or polyethylene (PE) material. This product may include specific features personalized to patients' retention plan. These features, selected by their doctor, may include but are not limited to, maintaining teeth position, maintaining jaw position, assisting with elastic wear, or helping manage oral habits. This product is not sterile and should be rinsed before placement in the mouth.

3. Product Performance

When used as prescribed, the product covers the crowns of patients' teeth. By matching the shape of the retainer to the position of patients' teeth, it helps hold them in place and maintain the results achieved after orthodontic treatment.

4. Directions for Use

The personalized treatment plan is based on the prescription. The following instructions should be reviewed with the patient.

a. Wearing procedure

- 1) Clean retainer: A soft-bristled toothbrush is recommended to clean the retainer. Dip the toothbrush in water to wipe the inner and outer surfaces of the retainer. Then rinse with water.
- 2) Clean mouth: Brush teeth and rinse mouth with water.
- 3) Clean hands: Wash hands before handling retainer.
- 4) Insert retainer: Make sure to have the correct retainer. The "U" denotes upper teeth or "L" for lower teeth. Fit retainer on anterior teeth first, then posterior teeth.
Procedure:
 - First, the anterior area of retainer is pressed onto the anterior teeth area.
 - Then the posterior area of retainer is pressed onto the posterior teeth area.
 - After the retainer is fit in, various types of occlusion can be conducted to perfectly fit the retainer with teeth.If part of retainer or retainer is not perfectly fit, the patient can use this part of teeth to bite medical cotton roll or chewies, and help the retainer in good fit.

b. Removing procedure

Principle: Remove the retainer from posterior teeth to anterior teeth, from lingual to buccal region, and alternately on both sides.

Procedure:

- 1) Maxillary retainer: The maxillary retainer can be removed from one side to the other. First, clamp the lingual side of one side's posterior retainer with fingers. Then pull the lingual edge of posterior retainer downward and buccally simultaneously. Both lingual and buccal side of retainer are

- pulled alternately until the retainer is detached from the teeth.
- 2) Mandible retainer: The mandible retainer can be removed either from one side to the other or both sides simultaneously. First, clamp the buccal side of one side's posterior retainer with fingers. Then pull the buccal edge of posterior retainer upward and buccally simultaneously. Both lingual and buccal side of retainer are pulled alternately until the retainer is detached from the teeth.
- c. Method of preservation: The retainer should be removed during eating, brushing teeth and can be removed for short social occasions. The retainer should be kept in a clean case after cleaning. It needs to be cleaned again before use.
- d. Wearing time
 - 1) The retainer should be worn as directed by the doctor.
 - 2) Before moving to a new retainer, present retainer should fully fit over the teeth, be easily removed and no orthodontic force should be felt. It is recommended to change into new retainer at night prior to bedtime.
 - 3) If the patient is unable to wear the retainer as prescriptive information, the patient should contact their doctor as soon as possible. To adjust the prescriptive wearing time, the doctor's consent must be obtained to continue the treatment. The treatment plan shall be followed, or the treatment shall be deemed as self-termination.
- e. Request for revisit
 - 1) The patient shall follow a regular revisit schedule according to the doctor's request to investigate the function and progress of the treatment.
 - 2) The current retainer should be worn at the time of revisit.

5. Safety Considerations

- a. The retainer must be used under the guidance of the attending doctor.
- b. The retainer is made for a particular patient and can only be used for that particular patient.
- c. If there is discomfort when wearing the retainer, patient should not adjust the retainer, but contact the doctor as soon as possible.
- d. After receiving the retainer, the patient shall wear it immediately to prevent unintended tooth movement.
- e. If the patient is unable to wear the retainer for the recommended time, does not comply with the doctor's instructions for wear, cannot complete revisits on time, or teeth erupt during treatment, it may affect the treatment maintenance after the completion of orthodontic treatment.
- f. Some medical conditions and the use of drugs may affect the treatment.
- g. Unconventional shapes, sprouting or missing teeth may affect the wearing of retainer or may affect the treatment maintenance.
- h. In cases of traction accessories, the patients are required to cooperate with the traction according to the doctor's instructions.
- i. The functional accessories may fall off due to improper intraoral habits, malocclusion, improper usage, etc., or in rare occasions may be swallowed, if this happens, patient should consult the doctor.
- j. In rare instances, problems in the temporo-mandibular or jaw joint (temporomandibular joint disorders (TMD)) may result in joint pain, headaches, or ear problems. Consult the doctor if these conditions occur during use.
- k. Patient should avoid improper intraoral habits, improper usage, etc. If the retainer is fractured or damaged, or the surface of the retainer is peeled, patient should contact doctor immediately.
- l. If the retainer cannot be engaged, contact the doctor immediately.
- m. If patient gets gum injury, oral ulcers, or abnormal loosening of teeth during use, patient should contact the doctor as soon as possible. If necessary, stop using the product and start again once the tissues heal.
- n. Failure to maintain good oral hygiene or consuming sugary or acidic foods and drinks while wearing the retainer may result in tooth decay (caries), gum disease (periodontal disease), staining, or enamel demineralization (white spots).
- o. When the retainer is damaged, deformed or lost, patient should contact the doctor and take prompt remedial measures.

INSTRUCTIONS FOR USE

6. Precautions

Retainer should be removed for eating, drinking and oral care.

After eating or drinking sugary drinks, patient should take oral cleaning measures (such as gargling with water, brushing their teeth, etc.) before wearing the retainer.

Patient should contact the doctor as soon as possible for the excessive irritation to the oral mucosa caused by wearing the retainer, and no self-adjustment by patient shall be allowed.

7. Indications for Use

Retainer is used to maintain orthodontic treatment results. By holding the teeth in their intended position.

8. Contraindications

Retainer is contraindicated for the patients with:

Serious progressive periodontal disease, or known TPU/PE allergies.

9. Residual risk

The residual risk has been judged as being acceptable. User-relevant safety information is provided in the Instructions for Use document of the medical device.

10. Maintenance

After each use, a soft tooth brush to clean the inside and outside surfaces of the retainer, and then put the retainer into a clean case for storage.

11. Storage Conditions

The product should be stored in a non-corrosive gas, cool and clean environment.

12. Production Date

See product packaging label for details.

13. Product Wearing Cycle

Wear this product for the duration prescribed by the doctor. It is recommended to wear each set for no more than 3 months.

14. Term of Use

Under proper storage conditions, the retainer has a shelf life of 3 years after production.

15. Label Description

The "U" in the label represents the upper retainer and is worn on the upper jaw; "L" means the lower retainer and is worn on the lower jaw.

16. Disposal and Treatment

The retainer is eco-friendly and does not contain environmentally harmful substances. Retainer should be disposed of in compliance with national and local laws.

*Notice: If any serious incident occurs in relation to the device, please report to the manufacturer and the competent authority of the Member State in which you are established.

17. Symbols glossary

Symbol	Descriptions	Symbol	Descriptions
	Batch Code		Use-by date
	Caution		Consult instructions for use or consult electronic instructions for use
	Medical Device		Manufacture Date
	Manufacturer Name and Address	Rx only	Prescription Use Only
	Model Number		Non sterile
	Patient Number		Unique device identifier
	Importer		Single patient multiple use

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Wuxi EA Bio-Tech Co., Ltd.